


BiteWing-Ease™ X-Ray Comfort Barriers

Large | White **REF** BWEL / Small | Blue **REF** BWES


INSTRUCTIONS FOR USE

- **Large | White:** Adjustable to fit all size 2 sensors
- **Small | Blue:** Adjustable to fit all size 0-1 sensors


Photos used in demonstration are BiteWing-Ease™ Comfort Barriers Small, Blue.



1
Insert sensor image-side up in comfort barrier.



2
Remove "Peel & Wrap" liner and fold to center to create bite tab.



3
Take images and discard comfort barrier after use.



Crosstex International, Inc., a Cantel Medical Company
6789 W. Henrietta Road, Rush, NY 14543 USA
585.359.0130 | 631.582.6777 | www.crosstex.com



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Not made with natural rubber latex.

Made in China
BWE-IFU-EN / REV B - 04/2021

All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.
©2021 Crosstex International, Inc. All rights reserved.