Steam Distiller

Purification System, Model DS1000

Pure Clean Water, Made Simple...

Simple one touch operation.

You will love using our model DS1000 Steam Distiller.

Fill the stainless-steel boiler with tap water and press the button. Shortly after you will have the purest possible water.

Features:

- Produces 1 gallon of 100% steam distilled water in 4 hours (Up to 6 gallons per day)
- Reduces 99% of total dissolved solids (TDS)
- Dripless stop-and-serve feature





PN# DS1000

How It Works:

The Tuttnauer Steam Distiller Model DS1000 system produces purified water by steam distillation. As the water temperature rises to 212°F (100°C), it kills bacteria, cysts, and viruses that may be present.

- 1. As the boiled water
- 2. Converts to steam
- 3. It rises and leaves behind dead microbes, and other contaminants and low boiling light gases are discharged through the gaseous vent
- 4. Steam is then cooled in the stainless-steel condenser turning it into high purity water.
- 5. Distilled water percolates through a final polishing coconut shell carbon filter
- 6. It is then collected in the Tritan™ carafe which is BPA free.



More Features:

- Automatic shutoff after each 1 gallon cycle
- 800 watts of power, 120 volt 60 hZ, fits standard outlet
- Effectively removes arsenic, nitrates, sodium, chlorine, bacteria and viruses.
- 1-gallon (4 liter) BPA Free carafe included
- Item Weight 11.43 pounds
- Dimensions: 15.25" H x 15" W x 9.75" D



The seamless, stainless steel boiler has a concealed heating element that won't become coated with scale and mineral deposits. A simple twist locks the 7 inch stainless steel cover firmly into place.

Accessory Items:



6-Pack Certified Carbon Filter Bags PN# 8606



Extra Tritan™ BPA-FREE Carafe PN# 3250



Kleenwise™ Distiller Cleaner and De-Scaler PN# KSD



Learn from our experts.

Visit our blog: tuttnauer.com/blog 🛭 🗓 🕞 🕞







Tuttnauer USA Co. Ltd.

25 Power Drive, Hauppauge, NY 11788

Tel: +800 624 5836, +631 737 4850, Fax: +631 737 0720 E-mail: info@tuttnauerUSA.com, www.tuttnauerUSA.com