

Monarch PSP Cleaning Wipes



Now
with softer
towelette
fabric

Extend the life of your PSPs and maintain optimal image quality with our specially formulated cleaning wipes.

- Safely remove dust, dirt and smudges from the imaging surface to maintain optimal image quality.
- Extra soft towelette fabric will not scratch or damage the plate.
- Convenient single-use wipes for fast, easy cleanup.
- Doesn't leave any residue. Evaporates completely.
- For all size PSPs and plate protectors.
- Economical: 50 wipes per box.

Cleaning Phosphor Storage Plates

Handle PSPs with care and keep them clean to achieve best image results.



1. Gently wipe the PSP Cleaning Wipe over the dry plate surface. Wipe back and forth and then in a circular motion.
2. Allow the plate surface to air dry. Make sure that the PSP is completely dry before re-using and inserting into the barrier envelope.

P/N B8910

 **Monarch**

by

AIR
TECHNIQUES equipped for life®

Visit airtechniques.com/product/psp-cleaning-wipes/ for more information.

Or contact our Monarch Support Team at 800-247-8324 ext. 5513

Order online at airtechniques.com/shop-list

PN 9922-558 Rev. A

Surface Disinfection

Instrument Cleaning

Skin & Hand Protection

Equipment Maintenance