

# Distinct


Aurelia® Distinct™ Honeycomb® Textured  
Powder Free Latex Examination Gloves



**Unique Honey Comb Texture  
Superior Tactile**



- Patented Honeycomb texture for extraordinary gripping
- Polymer-coating enhances donning and keeps hands smooth
- Excellent tactile sensitivity



**AURELIA®**  
Sign of Comfort™

[www.aureliagloves.com](http://www.aureliagloves.com) | 1-877-AURELIA (1-877-287-3542)

# Distinct

## Product Description

Aurelia® Distinct™ Honeycomb® Textured Powder Free Latex Examination Gloves are made from 100% natural Malaysian rubber, ambidextrous, creamy white color, Honeycomb® Textured and non-sterile. Aurelia gloves conform and are tested to the highest international standards. For more information visit [www.aureliagloves.com](http://www.aureliagloves.com)



## Aurelia® Distinct™ Specifications

### TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	970
Small	31.1	970
Medium	31.0	970
Large	31.2	960
Extra Large	31.3	970
ASTM 3578:05	18 MPa minimum	650% minimum

### AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	28.0	900
Small	28.0	920
Medium	27.5	910
Large	27.3	890
Extra Large	27.9	900
ASTM 3578:05	14 MPa minimum	500% minimum

## DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
<b>Extra Small</b> Aurelia Distinct ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>Small</b> Aurelia Distinct ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>Medium</b> Aurelia Distinct ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>Large</b> Aurelia Distinct ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>Extra Large</b> Aurelia Distinct ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum

### ITEM NUMBER

Extra Small	Small	Medium	Large	Extra Large
29225	29226	29227	29228	29229

## Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.



Authorized Aurelia Distributor:

CAUTION STATEMENT: This product contains natural rubber latex which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

AURELIA Div. Supermax Inc.

1899 Sequoia Drive, Aurora, IL 60506 | Toll Free: 1-877-AURELIA | Tel: 1-630-898-8886 | Fax: 1-630-898-8855  
Web: [www.aureliagloves.com](http://www.aureliagloves.com) | Email: [info@aueliagloves.com](mailto:info@aueliagloves.com)