

NEW!



Dryz[®] blu

Gingival Hemostatic Retraction Paste

Create a dry working environment every time!

- Features a bright aqua blue color allowing easy detection of where material is placed.
- Provides fast, effective hemostasis due to its aluminum chloride content.
- Stiff consistency for more displacement of sulcular tissue.
- Delivered in single dose capsules for more precise placement and asepsis.

New Dryz Blu features a bright, contrasting aqua blue color making the material more visible against tooth structure. The slight stiffness of the paste provides better displacement of tissue when the material is pressed subgingivally (with either a compression cap or cotton roll).

The material rinses easily and because of the aqua blue color, any residual is quickly spotted and removed simply with a water/air spray.

SKU #	Product
S186	NEW Dryz[®] Blu Intro Kit Includes (30) unit dose capsules (0.25 g/capsule) & Dispensing Gun.
S187	NEW Dryz[®] Blu Unit Dose Capsules Includes (30) unit dose capsules (0.24 g/capsule).
S190	NEW Dryz[®] Blu Syringes Includes (7) 0.5ml syringes and (15) applicator tips.
S191	NEW Dryz[®] Blu Syringes Value Pack Includes (25) 0.5ml syringes and (55) applicator tips.



parkell

Parkell.com
1-800-243-7446